

AH401B Picoammeter

For applications with multi-channel simultaneous acquisition



4-channel, 20-bit resolution, charge integration based, low noise digital picoammeter

HIGHLIGHTS

FEATURE

BENEFIT

Particularly suited for applications where multi-channel simultaneous acquisition is required, i.e. 4-quadrant photodiodes for beam displacement measurements

Performs current measurement from 50 pA (with a resolution of 50 aA) up to 1.8 μ A (resolution of 1.8 pA), with integration time ranging from 1 ms up to 1 s

Extremely low current measurements

Housed in a light and extremely compact box

Can be placed close to the signal sources in order to reduce cable lengths and minimize possible noise pick-up

Each input channel has two integrator stages

Current-to-voltage conversion can be performed continuously during the ADC conversion avoiding any dead time in the data output

Low temperature drifts, good linearity and very low noise

High-precision current measurements

Easily controlled via communication interface

Integration time, range, data format, type of acquisition, baud rate and many other parameters can be instantly set and checked

Modular communication capability

Allows the user to liberally select the type of communication interface, allowing control of the instrument with different types of programming languages and/or operating systems

APPLICATIONS

- Ultra-low current measurements
- Beam position monitoring
- Diamond detectors readout
- Ion chamber readout
- Laser diagnostics

HOW DOES IT WORK?

The compact user-friendly AH401B Picoammeter covers the whole acquisition, counting and digitization chain. To a great extent this simplifies and streamlines the read out at high speed and with extremely low noise X-ray detectors.

It is composed by a particular charge-integration input stage for low-current sensing combined with a 20-bit sigma-delta ADC converter integrating a noise reduction digital filter.

The acquisition of samples from the AH401B may be performed using either the “continuous” or “on demand” transmission modes:


- “continuous” mode: data are continuously sampled and transmitted, without external intervention, to the host device, allowing for real time data acquisition;
- “on demand” mode: data are sampled and transmitted only on a specific remote command request.

The external TRIGGER/GATE input signal is available to synchronize the acquisition of the picoammeter with external events (i.e. laser triggering). Furthermore, digital samples may be transferred using either the ASCII format or RAW binary data format for fast data transmission.

SPECIFICATIONS

| | |
|--|--|
| Input channels | 4 |
| Input connectors type | BNC |
| Effective current measuring range | From 50 pA to 1.8 μ A |
| Resolution | 20 bit |
| Data transfer | Up to 1 k samples/sec |
| Integration time | from 1 m sec to 1 sec |
| Polarity | Positive |
| Communication modules (AH401B) | Ethernet TCP-IP/UDP, USB 2.0, RS-232/422/485 |
| I/O Signal (AH401B) | CONV output - TRIGGER/GATE input |
| Supply voltage | From 9 V to 15 V |
| Supply current | From 100 to 350 mA depending on comm. module |
| Dimensions | 140 x 110 x 28 mm |
| Weight | 420 g |

AH401B BLOCK DIAGRAM



DELIVERABLES

- AH401B Picoammeter
- Power supply integration PS-2209
- Preinstalled Ethernet communication module
Other compatible modules are: RS232, RS422/485, USB and Ethernet (TCP/IP and UDP)
- Oscilloscope LabView Software

Contact us!

Industrial Liaison Office

Elettra - Sincrotrone Trieste S.C.p.A.
S.S. 14 - km. 163.5 in Area Science Park, 34149 Basovizza - Trieste, Italy
Tel. +39 040 3758303 - Fax +39 040 3758623
ilo@elettra.eu - <http://ilo.elettra.eu>



Elettra Sincrotrone Trieste